



MStim Arth LGT-235

Portable Electro-Stimulation Therapy Device

Technical Specifications		
	CE Version	FDA Version
Mode	TENS	TENS, MCR
Pulse Duration	50-400 μ s	50-500 μ s
Pulse Frequency	1-120Hz	1-120Hz
Intensity	0-60mA	TENS, NMES: 0-100mA MCR: 0-1000 μ A
Time	1-60min	1-60min
Battery Life	More than 4 hours	More than 4 hours

Longest

Guangzhou Longest Science & Technology Co., Ltd.
Address: 5&6F, Building B4, No.11, Kaiyuan Avenue, Science City,
Guangzhou Hi-tech Industrial Development Zone, 510530 Guangzhou, China
Tel: 86-20-66353999 Fax: 86-20-66353920
E-mail: export@longest.cn Website: www.longest.cn



Health Life Longest Care

Longest

MstimArth LGT-235

MstimArth LGT-235, works with low frequency electro stimulation (MCR and TENS). It is latest technology pain management stimulator for treatment inflammation and pain with knee arthritis.

Applications

Joint cold



Joint sprain



Tired



Active injury




Features


Both Knees Treated Simultaneously

Wireless Bluetooth Connection

Quick connection, easy and convenient



Knee Brace with Comfortable Conductive Fabric



2 Modes

TENS -- instant pain management at short time
MCR -- long-term use to repair damaged tissue