

Technical Specification

Chambers	3
Compression mode	4
Time range	1min-24hours or continuous
Pressure range	40-150mmHg
Hold time	0 - 10s (Step by 1s)
Interval time	10 - 60s (Step by 1s)
Interface	5" touch screen
Rechargeable battery	Yes (Lithium battery : 11.1V 9000mAh)
Adapter	Yes
Input voltage	AC 100 – 240 V 50 – 60 Hz
Dimension	250 (L) x 184 (W) x 150 (H) mm
Weight	2.5 kg

Certification



National Medical Products Administration



CE (MDD)



ISO13485 / MDSAP



INMETRO

Longest

Guangzhou Longest Science & Technology Co., Ltd.

Address: 5&6F, Building B4, No.11, Kaiyuan Avenue, Science City,
Guangzhou Hi-tech Industrial Development Zone, 510530 Guangzhou, China

Tel: 86-20-66353999 Fax: 86-20-66353920

E-mail: export@longest.cn Website: www.longest.cn



LGT-2200DVT

Compression Therapy Device



Health Life Longest Care


Longest




LGT-2200DVT is a special Pressotherapy system that is designed to decrease the risk of DVT(Deep Venous Thrombosis) and Pulmonary Embolism(PE)before, during and after surgery

4 compression modes with reusable & disposable sleeves meet all your need

⊕ Compression Modes



⊙ For Disposable Sleeve



*Notice: For one - leg treatment, please use single - head hose.


⊙ For Reusable Sleeve



*Notice: For two - leg treatment, please use double - head hose.

⊕ Introduction of deep vein thrombosis (DVT)

⊙ Reason



Deep vein thrombosis is one of the most common complications of bed rest, prevention is urgently needed

⊕ Features

- ⊙ 4 Compression Modes
- ⊙ 5" Touch Screen
- ⊙ Reusable & Disposable Sleeves
- ⊙ Gradient Pressure Setting
- ⊙ Intermittent Pneumatic Compression
- ⊙ Pressure of Every Chamber Adjustable
- ⊙ Rechargeable Battery
- ⊙ Powered by AC or DC
- ⊙ Automatic Memory for Last Setting

