

Technical Specification

Chambers	4
Compression mode	6
Time range	1min - 99min
Pressure range	30 - 150mmHg
Rechargeable battery	Yes (Battery life: more than 6 hours)
Adapter	Yes (15V 6A)
Input voltage	AC 100 - 240V 50 - 60Hz
Dimension	270(L)X148(W)X128(H)mm
Weight	1.65KG (without battery)

Certification

510(k) | K191862



CE(MDD)



Health Canada



ISO13485 / MDSAP



Therapeutic Goods Administration

Longest

Guangzhou Longest Science & Technology Co., Ltd.
Address: 5&6F, Building B4, No.11, Kaiyuan Avenue, Science City,
Guangzhou Hi-tech Industrial Development Zone, 510530 Guangzhou, China
Tel: 86-20-66353999 Fax: 86-20-66353920
E-mail: export@longest.cn Website: www.longest.cn



LGT-2200SP

FAST RECOVERY WITH COMPRESSION PUMP



Health Life Longest Care

Longest


SPEED UP THE RECOVERY



Quality Overlapping Sleeve




Quiet Operation



World compatible adapter




Multi-flush modes for best recovery



Light control unit



Wide recovery Flush range



Digital control and quick operation



Multi-options for sleeves



Up to 6500mAh rechargeable battery



PRODUCT OVERVIEW

LGT-2200SP is great for the user looking to get recovery while being mindful to budget. It allows you to apply muscle relax and recovery to one or two parts of your body at a time. Includes: 1 Pre-Charged Control Unit, High Quality Sleeves (depends on orders), 1 Set of Connection Hose, 1 World Compatible Adapter.

INDICATIONS FOR USE

The LGT-2200SP is an air compression systems indicated to temporarily relieve muscle tension and to increase blood circulation to get fast recovery to the treated areas.

RISKS OF THE SYSTEMS

The risks and benefits of using the LGT-2200SP are the similar as having manual massage. If the pressure feels uncomfortable, you can reduce the intensity or stop the session.

